

Monopolar Electrodes

Enhanced product performance in healthcare means best medical performance

We provide clinicians with the highest quality and safest German-designed products that equate to improved and efficient delivery of healthcare.

Surgica offers an excessive range of high-quality Monopolar electrodes, perfect for your cutting, blending, desiccations and fulguration needs.

Our extensive range covers loop electrodes (suitable for LEEP), ball electrodes (crafted for the ablation/coagulation of soft tissue), a variety of needle electrodes (for every diagnosis purpose you need) and blade electrodes.

They team perfectly with our range of smoke evacuation pencils, standard disposable and reusable electrosurgical pencils, which are available in Chapter Four.

Browse below for our selection of reusables, stainless steel disposables and non-stick disposables.

REUSABLE RANGE	DISPOSABLE RANGE
Knife	Blade
Blade	Modified Blade
Modified Blade	Precision Blade
Needle	Ultra Fine Blade
Modified Needle	Needle
Tungsten Needle	Modified Needle
Ball	Blunt Needle
Modified Ball	Tungsten Needle
Lletz Loop	Ball
	Bayonet
	Lletz Loop
	Conization

nonstick

Our non-stick range are insulated, providing excellent thermal qualities with a reduced need to clean them during coagulation.

The non-stick tips give **increased accuracy** and **speed of dissection** with minimal tissue adherence or charring.

01 Reusable Monopolar Electrodes

Range covers: Knife, Blade, Modified Blade, Needle, Modified Needle, Tungsten Needle, Ball, Modified Ball, Lletz Loop.

Available for blade, ball, loop and needle electrodes

Value for money and cost efficiency

Can be used on lower power settings (controls thermal necrosis)

Comes in an array of sizes and styles

Surgica's Monopolar electrodes are the must-have tools for any practitioner involved in a heavier surgery. Designed in Germany, our Stainless-Steel tools are compatible with all speciality pencils.

Autoclavable, our tools can survive sterilisation at temperatures up to 138° without damage or peeling. With a faster start-up time and a lower

core temperature, our electrodes are designed to meet the surgeon's quick expectations – time wasted is a life lost.

Our Surgica range covers a variety of electrode possibilities, making it easier for you to find the best tool for the job.

Measuring instructions



Made from high quality stainless steel

Available in two shaft diameters (i.e. 2.38mm Ø and 4mm Ø)

Autoclavable at temperatures up to 138°

Monopolar Electrodes

Reusable

Knife Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP WIDTH	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP WIDTH
03-001	03-151	40mm	17mm	1.5mm	03-004	03-154	40mm	24mm	3.4mm
03-002	03-152	40mm	17mm	1.5mm	03-005	03-155	120mm	12mm	2.4mm
03-003	03-153	40mm	18mm	3.4mm	03-006	03-156	115mm	12mm	2.4mm

Blade Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP WIDTH	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP WIDTH
03-007	03-157	45mm	19mm	2.3mm	03-011	03-161	120mm	19mm	2.3mm
03-008	03-158	45mm	19mm	2.3mm	03-012	03-162	115mm	19mm	2.3mm
03-009	03-159	45mm	24mm	3mm	03-013	03-163	120mm	12mm	2.4mm
03-010	03-160	45mm	24mm	3mm	03-014	03-164	180mm	19mm	2.3mm

Modified Blade Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP WIDTH	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP WIDTH
03-015	03-165	45mm	19mm	2.3mm	03-017	03-167	120mm	19mm	2.3mm
03-016	03-166	45mm	19mm	2.3mm	03-018	03-168	180mm	19mm	2.3mm

Needle Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP DIAMETER
03-019	03-169	35mm	17mm	0.8mm	03-023	03-173	120mm	15mm	0.8mm
03-020	03-170	35mm	20mm	0.8mm	03-024	03-174	120mm	10mm	0.8mm
03-021	03-171	40mm	22mm	0.8mm	03-025	03-175	180mm	15mm	0.8mm
03-022	03-172	65mm	29mm	0.8mm					

Modified Needle Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP DIAMETER
03-026	03-176	40mm	22mm	0.8mm	03-027	03-177	115mm	15mm	0.8mm

Monopolar Electrodes

Reusable

Tungsten Needle Electrode

5 95 °C 138 °C 1kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP LENGTH	TIP DIAMETER
03-028	03-178	40mm	3mm	0.5mm	03-030	03-180	55mm	3mm	0.5mm
03-029 (angled)	03-179 (angled)	35mm	3mm	0.5mm	03-031 (angled)	03-181 (angled)	55mm	3mm	0.5mm

Ball Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP DIAMETER
03-032	03-182	40mm	2mm	03-037	03-187	35mm	4mm
03-033	03-183	40mm	3mm	03-038	03-188	110mm	3mm
03-034	03-184	40mm	4mm	03-039	03-189	110mm	4mm
03-035	03-185	40mm	5mm	03-040	03-190	110mm	5mm
03-036	03-186	35mm	2mm				

Modified Ball Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	TIP DIAMETER
03-041	03-191	115mm	2mm	03-042	03-192	115mm	4mm

Wire Snare Electrode

5 95 °C 138 °C 4.3kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	LOOP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	LOOP DIAMETER
03-043	03-193	40mm	6mm	03-046	03-196	105mm	5mm
03-044	03-194	40mm	12mm	03-047	03-197	105mm	10mm
03-045	03-195	40mm	16mm				

Square Loop Electrode

5 95 °C 138 °C 1.5kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	LOOP WIDTH	LOOP HEIGHT	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	LOOP WIDTH	LOOP HEIGHT
03-048	03-198	125mm	5mm	5mm	03-050	03-200	135mm	10mm	12mm
03-049	03-199	125mm	6mm	15mm					

Round Loop Electrode

5 95 °C 138 °C 1.5kVp



Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	LOOP DIAMETER	Shaft Ø 2.38mm	Shaft Ø 4.0mm	WORKING LENGTH	LOOP DIAMETER
03-051	03-201	130mm	10mm	03-055	03-205	140mm	20mm
03-052	03-202	135mm	15mm	03-056	03-206	145mm	25mm
03-053	03-203	135mm	15mm	03-057	03-207	150mm	30mm
03-054	03-204	140mm	20mm				

02 Disposable Monopolar Electrodes

Range covers: Blade, Modified Blade, Precision Blade, Ultra Fine Blade, Needle, Modified Needle, Blunt Needle, Modified Blunt Needle, Tungsten Needle, Ball, Bayonet, Lletz Loop, Conization.

Easily disposable for maximum hygiene where it's needed most


Value for money and cost efficiency

Can be used on lower power settings (controls thermal necrosis)

The disposable version of our high-quality monopolar electrodes, Surgica's brand offers the highest quality available for our tool's one-time use. Our stainless-steel tools are made for the tougher surgery, lasting through every dissection, coagulation and fulguration. Pre-sterilized, our tools come packaged to avoid damage during transit, saving valuable time.


Surgica's disposable monopolar electrodes offer the same quality as our reusable electrodes, including our in-built working mechanism and ergonomic handle, minimizing fatigue and working with the practitioner to maximize life-saving performance.

Measuring instructions



Designed in Germany

Created with high quality stainless steel



Two types of electrodes (i.e. with and without hexlock)

Comes in an array of sizes and styles

Monopolar Electrodes

Disposable

Blade Electrode

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-301	03-451	63.5mm (2 1/2")	44.3mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
								
03-302	03-452	70mm (2 3/4")	50.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
								
03-303	03-453	102mm (4")	82.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
								
03-304	03-454	165mm (6 1/2")	145.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
								
03-305	03-455	205mm (8")	185.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
								

Modified Blade Electrode

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-306	03-456	63.5mm (2 1/2")	44.3mm	4mm	2.55mm	12	Sterile EO	4.3kVp
								
03-307	03-457	70mm (2 3/4")	50.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
								
03-308	03-458	102mm (4")	82.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
								

continued...

Monopolar Electrodes

Disposable



Modified Blade Electrode

modified **BLEDA**

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-309	03-459	165mm (6 1/2")	145.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
								
03-310	03-460	205mm (8")	185.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
								

Blade Electrode

precision

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-311	03-461	70mm (2 3/4")	50.8mm	13.5mm	1.30mm	12	Sterile EO	4.3kVp
								
03-313	03-463	102mm (4")	82.8mm	13.5mm	1.30mm	12	Sterile EO	4.3kVp
								

Modified Blade Electrode

precision

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-312	03-462	70mm (2 3/4")	50.8mm	4mm	1.30mm	12	Sterile EO	4.3kVp
								
03-314	03-464	102mm (4")	82.8mm	4mm	1.30mm	12	Sterile EO	4.3kVp
								

Monopolar Electrodes

Disposable

Blade Electrode

ultra*line*

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-315	03-465	70mm (2 3/4")	50.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-317	03-467	102mm (4")	82.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-319	03-469	165mm (6 1/2")	145.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp

Modified Blade Electrode

ultra*line*

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-316	03-466	70mm (2 3/4")	50.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
03-318	03-468	102mm (4")	82.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
03-320	03-470	165mm (6 1/2")	145.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp

Monopolar Electrodes

Disposable

Needle Electrode

needle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-321	03-471	70mm (2 3/4")	50.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp
03-322	03-472	102mm (4")	82.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp
03-323	03-473	152mm (6")	132.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp
03-324	03-474	205mm (8")	185.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp

Modified Needle Electrode

modified*needle*

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-325	03-475	70mm (2 3/4")	50.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp
03-326	03-476	102mm (4")	82.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp
03-327	03-477	152mm (6")	132.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp
03-328	03-478	205mm (8")	185.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp

Monopolar Electrodes

Disposable

Blunt Needle Electrode

blunt needle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER				
03-329	03-479	70mm (2 3/4")	50.8mm	17.5mm	0.9mm	12	Sterile EO	4.3kVp	
03-330	03-480	70mm (2 3/4")	50.8mm	4mm	0.9mm	12	Sterile EO	4.3kVp	
Modified Blunt Needle Electrode									
03-331	03-481	102mm (4")	82.8mm	17.5mm	0.9mm	12	Sterile EO	4.3kVp	

Needle Electrode

tungsten needle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER				
03-332	03-482	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp	
03-333	03-483	63.5mm (2 1/2")	44.3mm	4mm	0.5mm	12	Sterile EO	1kVp	
03-334	03-484	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp	
45 Degree									
03-335	03-485	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp	
45 Degree									
03-336	03-486	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp	
90 Degree									

continued...

surgica

Monopolar Electrodes

Disposable

Needle Electrode

tungsten needle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-337	03-487	165mm (6 1/2")	145.8mm	4mm	0.5mm	12	Sterile EO	1kVp

Ball Electrode

ball

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP DIAMETER			
03-338	03-488	50mm (2")	30.8mm	4.25mm	12	Sterile EO	4.3kVp
03-339	03-489	70mm (2 3/4")	50.8mm	4.25mm	12	Sterile EO	4.3kVp
03-340	03-490	127mm (5")	107.8mm	4.25mm	12	Sterile EO	4.3kVp

Bayonet Electrode

bayonet


HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH/DIAMETER				
03-341	03-491	146mm (5 3/4")	126.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp	
03-342	03-492	146mm (5 3/4")	126.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp	
Modified Needle									
03-343	03-493	146mm (5 3/4")	126.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp	
Modified Blade									

surgica

SURQO LOOP

Disposable Monopolar Electrodes

Minimal loop flex for maximum performance



Designed for easy use and excellent cutting characteristics

Ideal for obtaining precise specimens for examination

Loop made of high-quality tungsten wire


Stable construction design, minimal loop flex

Comes in pack of 12

Precise tissue excision

Highly cost effective and VFM

Easy fit to electro-surgical pencils



Stable construction design, minimal loop flex

Colour-coded for shape and size recognition

Minimum output power settings (4.3kVp)

Monopolar Electrodes

Disposable

Square Loop Electrode

square

CODE	SHAFT LENGTH	LOOP WIDTH	LOOP HEIGHT			
03-344	120mm (4 3/4")	10mm	4mm	12	Sterile EO	1.5kVp
03-345	120mm (4 3/4")	10mm	7mm	12	Sterile EO	1.5kVp
03-346	120mm (4 3/4")	10mm	10mm	12	Sterile EO	1.5kVp

Radius Loop Electrode

radius

CODE	SHAFT LENGTH	LOOP WIDTH	LOOP HEIGHT			
03-347	120mm (4 3/4")	10mm	5mm	12	Sterile EO	1.5kVp
03-348	120mm (4 3/4")	10mm	10mm	12	Sterile EO	1.5kVp
03-349	120mm (4 3/4")	15mm	6mm	12	Sterile EO	1.5kVp
03-350	120mm (4 3/4")	15mm	12mm	12	Sterile EO	1.5kVp
03-351	120mm (4 3/4")	20mm	8mm	12	Sterile EO	1.5kVp









continued...

Monopolar Electrodes

Disposable

Radius Loop Electrode

radius








CODE	SHAFT LENGTH	LOOP WIDTH	LOOP HEIGHT			
03-352	120mm (4 3/4")	20mm	10mm	12	Sterile EO	1.5kVp
						
03-353	120mm (4 3/4")	20mm	12mm	12	Sterile EO	1.5kVp
						
03-354	120mm (4 3/4")	20mm	15mm	12	Sterile EO	1.5kVp
						
03-355	120mm (4 3/4")	25mm	8mm	12	Sterile EO	1.5kVp
						
03-356	120mm (4 3/4")	25mm	22mm	12	Sterile EO	1.5kVp
						

Monopolar Electrodes

Disposable

Conization Electrode

conization

CODE	SHAFT LENGTH	CONE RADIUS	HEIGHT			
03-357	120mm (4 3/4")	9mm	13mm	12	Sterile EO	1.5kVp
						
03-358	120mm (4 3/4")	12mm	10mm	12	Sterile EO	1.5kVp
						
03-359	120mm (4 3/4")	11mm	18mm	12	Sterile EO	1.5kVp
						
03-360	120mm (4 3/4")	15mm	23mm	12	Sterile EO	1.5kVp
						

03 *nonstick* Disposable Monopolar Electrodes

Range covers: Blade, Modified Blade, Precision Blade, Ultra Fine Blade, Needle, Modified Needle, Blunt Needle, Modified Blunt Needle, Tungsten Needle, Ball, Bayonet.

All product tips are coated in Polytetrafluoroethylene (PTFE), reducing flaking at high temperatures, minimizing friction during incision and allowing for a ready-to-wipe surface.

The tips reduce eschar build-up, lowering energy setting and saving money to the practitioner. Our tips are evenly coated to ensure consistent coagulation.

surqo clean

Our range of disposable easy-cleaning electrode equipment


Available for bayonet, blade, ball and needle electrodes

Easy disposal for maximum hygiene where it's needed most

Less eschar build-up so less cleaning

Value for money and cost efficiency

Measuring instructions



GREEN PTFE (without hexlock) GREY PTFE (with hexlock)

PTFE coating for a fast, time-saving clean-up

Two types of electrodes (i.e. with and without hexlock)

Minimum output power settings (4.3kVp)

Monopolar Electrodes

Non Stick Disposable

Blade Electrode

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-601	03-751	63.5mm (2 1/2")	44.3mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-602	03-752	70mm (2 3/4")	50.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-603	03-753	102mm (4")	82.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-604	03-754	165mm (6 1/2")	145.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-605	03-755	205mm (8")	185.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp

Modified Blade Electrode

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-606	03-756	63.5mm (2 1/2")	44.3mm	4mm	2.55mm	12	Sterile EO	4.3kVp
03-607	03-757	70mm (2 3/4")	50.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
03-608	03-758	102mm (4")	82.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp

continued...

Monopolar Electrodes

Non Stick Disposable

Modified Blade Electrode

modified blade

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-609	03-759	165mm (6 1/2")	145.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp



03-610	03-760	205mm (8")	185.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
--------	--------	------------	---------	-----	--------	----	------------	--------



Blade Electrode

precision

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-611	03-761	70mm (2 3/4")	50.8mm	13.5mm	1.30mm	12	Sterile EO	4.3kVp



03-613	03-763	102mm (4")	82.8mm	13.5mm	1.30mm	12	Sterile EO	4.3kVp
--------	--------	------------	--------	--------	--------	----	------------	--------



Modified Blade Electrode

precision

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-612	03-762	70mm (2 3/4")	50.8mm	4mm	1.30mm	12	Sterile EO	4.3kVp



03-614	03-764	102mm (4")	82.8mm	4mm	1.30mm	12	Sterile EO	4.3kVp
--------	--------	------------	--------	-----	--------	----	------------	--------



Monopolar Electrodes

Non Stick Disposable

Needle Electrode

needle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-621	03-771	70mm (2 3/4")	50.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp



03-622	03-772	102mm (4")	82.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp
--------	--------	------------	--------	--------	-------	----	------------	--------



03-623	03-773	152mm (6")	132.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp
--------	--------	------------	---------	--------	-------	----	------------	--------



03-624	03-774	205mm (8")	185.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp
--------	--------	------------	---------	--------	-------	----	------------	--------



Modified Needle Electrode

modified needle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-625	03-775	70mm (2 3/4")	50.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp



03-626	03-776	102mm (4")	82.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp
--------	--------	------------	--------	-----	-------	----	------------	--------



03-627	03-777	152mm (6")	132.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp
--------	--------	------------	---------	-----	-------	----	------------	--------



03-628	03-778	205mm (8")	185.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp
--------	--------	------------	---------	-----	-------	----	------------	--------



Monopolar Electrodes

Non Stick Disposable

Blunt Needle Electrode

bluntneedle

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-629	03-779	70mm (2 3/4")	50.8mm	17.5mm	0.9mm	12	Sterile EO	4.3kVp
03-630	03-780	70mm (2 3/4")	50.8mm	4mm	0.9mm	12	Sterile EO	4.3kVp
03-631	03-781	102mm (4")	82.8mm	17.5mm	0.9mm	12	Sterile EO	4.3kVp

Bayonet Electrode

bayonet

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH/DIAMETER				
03-641	03-791	146mm (5 3/4")	126.8mm	17.5mm	0.8mm	12	Sterile EO	4.3kVp	
03-642	03-792	146mm (5 3/4")	126.8mm	4mm	0.8mm	12	Sterile EO	4.3kVp	
Modified Needle									
03-643	03-793	146mm (5 3/4")	126.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp	
Modified Blade									

Monopolar Electrodes

Non Stick Disposable

Ball Electrode

ball

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP DIAMETER				
03-638	03-788	50mm (2")	30.8mm	4.25mm	12	Sterile EO	4.3kVp	
03-639	03-789	70mm (2 3/4")	50.8mm	4.25mm	12	Sterile EO	4.3kVp	
03-640	03-790	127mm (5")	107.8mm	4.25mm	12	Sterile EO	4.3kVp	

surqo'ine

Our ultra-fine needle electrode range

Ultra-fine needle electrodes offer improved cutting and coagulation performance, while minimizing thermal spread. They offer surgeons the precision they need for fine dissection in a wide variety of procedures. Lower setting is just as effective for surgery, which helps control thermal necrosis and nearby tissue and nerve damage.

The PTFE coated needle reduces eschar build-up, allowing for quicker clean-up, whilst withstanding sterilisation temperatures with minimum damage.


Ultrafine tip decreases power settings, providing maximum power at energy saving costs

Controls thermal necrosis, minimising damage to external tissue

Reduces eschar build-up, providing an easier cleaning solution

Angular design available (45 or 90 degrees)

Measuring instructions



GREEN PTFE (without hexlock) GREY PTFE (with hexlock)

Coated with PTFE for non-stick coverage

Two types of electrodes
(i.e. with and without hexlock)

Sizes available from 50mm (2")
(tip length 4mm, diameter 0.5mm)

Monopolar Electrodes

Non Stick Disposable

Needle Electrode

ultra'ine

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP DIAMETER			
03-632	03-782	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp
03-633	03-783	63.5mm (2 1/2")	44.3mm	4mm	0.5mm	12	Sterile EO	1kVp
03-634	03-784	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp
45 Degree								
03-635	03-785	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp
45 Degree								
03-636	03-786	50mm (2")	30.8mm	4mm	0.5mm	12	Sterile EO	1kVp
90 Degree								
03-637	03-787	165mm (6 1/2")	145.8mm	4mm	0.5mm	12	Sterile EO	1kVp

surqoblade

Our non-stick, ultra-fine blade range

Surqoblade is the ultimate multi-functional scalpel blade, perfect for effective cutting, coagulation and dissection. Its efficient design means it works as efficiently as the best standard scalpel, whilst maintaining the same quality electrode work that the rest of the Surgica Range provides. The quality of dissection allows wound healing time equivalent of a standard scalpel, minimising tissue damage and providing a quick, clean cut.

Our stainless-steel blades effortlessly cut through tissue, with our PTFE coating providing non-stick protection, reducing eschar build-up and saving time with quick clean-up.

Our product provides a safer, more practical alternative to the standard scalpel, with all added benefit of being high-quality monopolar electrode.

A powerful electrode and practical alternative to the scalpel – sharp edges do the work for you

Precise design, built to handle incision, coagulation and dissection with maximum efficiency

Stainless steel blades to deal with even the most resilient tissue

PTFE coat, providing a non-stick resolution to clean-ups

Measuring instructions



GREEN PTFE (without hexlock) GREY PTFE (with hexlock)

Less eschar build-up,
higher functioning blade

Two types of electrodes
(i.e. with and without hexlock)

Sizes available from 70mm (2¾")
(tip length 4mm / 17.5mm, width 2.55mm)

Monopolar Electrodes

Non Stick Disposable

Blade Electrode

ultrafine

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-615	03-765	70mm (2 ¾")	50.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-617	03-767	102mm (4")	82.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp
03-619	03-769	165mm (6 ½")	145.8mm	17.5mm	2.55mm	12	Sterile EO	4.3kVp

Modified Blade Electrode

ultrafine

HEXLOCK	WO / HEXLOCK	OVERALL LENGTH	WORKING LENGTH	TIP LENGTH	TIP WIDTH			
03-616	03-766	70mm (2 ¾")	50.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
03-618	03-768	102mm (4")	82.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp
03-620	03-770	165mm (6 ½")	145.8mm	4mm	2.55mm	12	Sterile EO	4.3kVp