

Monopolar Forceps

Much like our bipolar range, **Surgica provides a wide range of monopolar forceps**, both reusable and disposable.

Our reusable range are made with the finest German stainless steel, with our disposable range made from the highest-grade stainless steel.

Both are designed to provide **control, safety and maximum surgical efficiency**, even at the lowest setting. Our extensive range are insulated, providing excellent thermal qualities with a reduced need to clean them during coagulation.

Our precise Monopolar tips give increased accuracy and speed of dissection with minimal tissue adherence or charring.

FORCEPS DESIGNS

- Adson
- McIndoe
- Gillies
- Potts-Smith
- Jansen Bayonet



surgica

Types of grip



SERRATED FLAT
The classic flat tips with serrations provide improved grip for grasping of tissues, excellent sealing of vessels, perfect for endoscope assisted procedures.



SERRATED CLEARVIEW
The short and fine design is combined with clearview adaption, providing an obstacle free view of the surgical area with serrated tips, ideal for precise grasping and manipulating of tissue during more precise surgeries..

Tip widths

All tips come in three styles.



1.0mm
fine



1.5mm
blunt thin



2.0mm
blunt thick



surgica

01 Reusable Monopolar Forceps

Designed in Germany from the finest stainless steel.

- Longer lasting lifespan than standard steel
- Quick coagulation in seconds, even at lowest settings
- Atraumatic design avoids unnecessary tissue disruption
- Made for tougher surgeries (harder cuts and coagulation)

Surgica’s standard monopolar forceps are made to help the surgeon during longer, tougher surgeries. Our top-of-the-range German inspired tools are professionally made with the finest, highest quality stainless steel. Its longer lasting lifespan isn’t affected by the more scrutinous surgeries, as these tools are built to last.

Our quick, easy to clean tools makes the post-surgery wipe down much easier, vastly reducing eschar build-up. Our tools can be sterilized n temperatures up to 138° without flaking or damage, allowing for a cleaner, safer antibacterial steam.

Much like our bipolar forceps, our Monopolar forceps feature fast start-up, reaching maximum temperature in seconds, with its core temperature being lower than the average Monopolar forcep.

Our slimline designed monopolar forceps are the perfect tool for the job. Its sleek design and in-built mechanisms aid the surgeon across all aspects of surgery, including tough cuts and coagulations. The atraumatic handle rests the surgeon’s hand, avoid hand strain, allowing the device to work for longer.



- Made from the finest German Stainless Steel
- Ergonomic handle minimises hand muscle fatigue
- Autoclavable at temperatures up to 138°


Monopolar Forceps

Reusable

Adson, 121mm (4 ¾”)

1 138 °C 2kVp


- 02-001 1.0mm, straight
- 02-002 1.5mm, straight



McIndoe, 152mm (6”)

1 138 °C 2kVp


- 02-007 1.5mm, straight
- 02-008 2.0mm, straight



Gillies, 152mm (6”)

1 138 °C 2kVp


- 02-013 1.5mm, straight
- 02-014 2.0mm, straight



Potts-Smith, 178mm (7”)

1 138 °C 2kVp


- 02-019 1.5mm, straight
- 02-020 2.0mm, straight



Potts-Smith, 215mm (8 ½”)

1 138 °C 2kVp

- 02-025 1.5mm, straight
- 02-026 2.0mm, straight



Monopolar Forceps

Reusable

Potts-Smith, 240mm (9 ½”)

1 138 °C 2kVp

02-031 1.5mm, straight | 02-032 2.0mm, straight



Jansen Bayonet, 165mm (6 ½”)

1 138 °C 2kVp

02-037 1.5mm, straight | 02-038 2.0mm, straight



Jansen Bayonet, 190mm (7 ½”)

1 138 °C 2kVp

02-043 1.5mm, straight | 02-044 2.0mm, straight | 02-049 1.5mm, angled up | 02-050 2.0mm, angled up



Jansen Bayonet, 215mm (8 ½”)

1 138 °C 2kVp

02-055 1.5mm, straight | 02-056 2.0mm, straight | 02-061 1.5mm, angled up | 02-062 2.0mm, angled up



02 Disposable Monopolar Forceps

Designed in Germany from the highest quality stainless steel.

Polished finish allows easier wipe

Quick coagulation at lower peak temperatures

Ergonomic handle and in-built mechanisms

Built for quick use on harder surgeries

Surgica’s standard monopolar disposable forceps are made for quick, one-time use in tougher surgeries.

Designed in Germany and **manufactured using the highest quality stainless steel**, our disposables share the same quality attributes as our standard reusables, including the extremely quick coagulation in seconds, the lower peak temperatures increasing the surgeons’ carbon footprint, and the better, stronger coagulation, built for tougher surgeries.

Our disposable monopolar forceps focus to make its one-time use the best a surgeon can have. With our in-built working mechanism and ergonomic

handle, it makes for comfortable, functional usage, minimizing fatigue. The lightweight design is made for quick use, whilst our tools come pre-sterile, perfect for situations that need fast response times and instant use.

The smooth stainless-steel finish makes **wiping fast and cleaning easier**. On top of all this, our monopolar disposables are easily disposed of, built for both European and American usage. You have the option of straight or bayonet, and our forceps comes with an integrated cable, making our tools essential for any surgery.

Measuring instructions



Plugs & connectors



Built for quick use on harder surgeries, deeper cuts and harder coagulations

Integrated cable with 8mm plug


Monopolar Forceps

Disposable

Adson, 121mm (4 ¾")

20 Sterile EO 2kVp


02-151 1.0mm, straight | 02-152 1.5mm, straight



McIndoe, 152mm (6")

20 Sterile EO 2kVp

02-157 1.5mm, straight | 02-158 2.0mm, straight



Gillies, 152mm (6")

20 Sterile EO 2kVp

02-163 1.5mm, straight | 02-164 2.0mm, straight



Potts-Smith, 178mm (7")

20 Sterile EO 2kVp


02-169 1.5mm, straight | 02-170 2.0mm, straight



Potts-Smith, 215mm (8 ½")

20 Sterile EO 2kVp

02-175 1.5mm, straight | 02-176 2.0mm, straight



Monopolar Forceps

Disposable

Potts-Smith, 240mm (9 ½")

20 Sterile EO 2kVp

02-181 1.5mm, straight | 02-182 2.0mm, straight



Jansen Bayonet, 165mm (6 ½")

20 Sterile EO 2kVp

02-187 1.5mm, straight | 02-188 2.0mm, straight



Jansen Bayonet, 190mm (7 ½")

20 Sterile EO 2kVp

02-193 1.5mm, straight | 02-194 2.0mm, straight | 02-199 1.5mm, angled up | 02-200 2.0mm, angled up



Jansen Bayonet, 215mm (8 ½")

20 Sterile EO 2kVp

02-205 1.5mm, straight | 02-206 2.0mm, straight | 02-211 1.5mm, angled up | 02-212 2.0mm, angled up

